Q1 series

FEATURES

Up to 40 mJ pulse energy, up to 0.5 W average output power

Up to 50 Hz variable pulse repetition rate

Air cooled (water-free)

7 – 10 ns pulse duration

Weight < **6 kg** incl. heatsink, controller & AC/DC adapter

> 2 G shot lifetime of pump diode

Built-in sync pulse generator for triggering of user equipment

Remote monitoring and control via built-in **Ethernet** interface

Optional attachable PC controlled motorized attenuator

Optional attachable pulse energy monitor

Optional fiber coupled output

Optional baseplate for mounting Q1 and H1 or Q1 and fiber coupler together

Optional stand-alone 2^{nd} , 3^{rd} , 4^{th} or 5^{th} harmonic generator

Stand-alone air-purging unit for long lifetime of UV optics

APPLICATIONS

Laser Induced Breakdown Spectroscopy (LIBS)

Laser ablation

Time-of-Flight Spectroscopy (TOFS)

Light Induced Fluorescence (LIF) Spectroscopy

Non-linear laser spectroscopy

Flash photolysis

Photoacoustics

Metrology



DIODE PUMPED AIR-COOLED Q-SWITCHED LASER



Q1 is compact, energy-efficient, diode pumped, air-cooled, Q-switched laser designed for wide range of applications that require high peak power pulses.

Laser can produce up to 40 mJ at 10 Hz pulse repetition rate or up to 10 mJ at 50 Hz pulse repetition rate. Low divergence of laser beam allows efficient conversion to harmonic wavelengths using optional H1 series harmonic generator module.

The laser can be configured to emit 1053 nm or 1064 nm wavelength from Nd:YLF or Nd:YAG laser crystals respectively. Due to thermal properties of Nd:YLF crystal, at 1053 nm the laser can operate from single shot to maximum pulse repetition rate without significant changes in beam divergence or profile.

Innovative laser design resulted in user-friendly turnkey system that requires little maintenance. There are no chillers or bulky power supplies that one needs fit under the table. All laser electronics are integrated into Q1 housing and the only external modules are lightweight controller box that provides laser control interfaces and mains adapter that provides 12 VDC powering. Total weight of the laser is less than 5 kg.

Laser is monitored and controlled trough Ethernet port via build-in web-server. Any computer or even cell phone with modern web-browser will be able to control Q1. API is provided for integration with user devices. Low jitter triggering pulses for user equipment are available with up to 300 μ s lead in internal triggering mode. In external triggering mode, laser pulsing can be externally triggered from delay generator.

Laser functionality can be further extended by wide selection of auxiliary equipment:

- > Up to fifth harmonic output can be produced from our stand-alone H1 series harmonic generators.
- > Lifetime of UV optics can be extended by stand-alone closed loop air-purging unit.
- > Pulse energy can be adjusted with our motorized attenuator for fundamental wavelength beam.
- > Pulse energy can be monitored by our attachable pulse energy monitor with analog and/or digital output.
- > Fiber coupled output is available by request. Please inquire for detailed specifications.

SPECIFICATIONS AT **10 Hz** PULSE REPETITION RATE ¹⁾

WODEL	Q1				
MODEL	-B10	-C10	-D10		
Wavelength	1064 or	1053 nm	1064 nm	1053 nm	
Pulse repetition rate 2)	single-shot to 10 Hz				
Pulse duration 3)		< 8	3 ns		
Pulse energy 4)					
1064 / 1053 nm	8 / 10 mJ	16 / 20 mJ	32 mJ	40 mJ	
532 / 526.5 nm	4 / 5 mJ	8 / 10 mJ	16 mJ	20 mJ	
355 / 351 nm	2.4 / 3 mJ	4.5 / 6 mJ	9 mJ	12 mJ	
266 / 263 nm	1.2 / 1.5 mJ	2.4 / 3 mJ	5 mJ	6 mJ	
213 / 211 nm	0.5 mJ	1 mJ	2 mJ	2.5 mJ	
Pulse to pulse energy stability 5)					
1064 / 1053 nm	< 0.5 % RMS				
532 / 526.5 nm	< 2.5 % RMS				
355 / 351 nm	< 3.5 % RMS				
266 / 263 nm	< 4 % RMS				
213 / 211 nm	< 5 % RMS				
Power drift 6)	± 3.0 %				
Beam profile	bell-shaped, >75 % fit to Gaussian				
Beam divergence 7)	< 1.5 mrad	< 1 mrad			
Beam pointing stability 8)	< 10 μrad				
Polarization		linear, horizontal			
Typical beam diameter 9)	1.7 mm	2.5 mm 3.5 mm		mm	
Jitter 10)		< 0.5 n	is RMS		
OPTIONAL ATTENUATOR 11)					
Transmission range	0.5 – 95 %				
DIMENSIONS					
Laser head (W×L×H)	110 × 231× 112 mm ³				
Controller unit (W×L×H)		108 × 191	× 59 mm ³		
Power adapter (W×L×H) 12)	$50 \times 125 \times 31 \text{ mm}^3$				
Harmonic generator (W×L×H) 13)	110 × 242 × 112 mm ³				
OPERATING REQUIREMENTS					
Cooling requirements	air cooled				
Ambient temperature	15 – 30 °C				
Relative humidity	10 – 80 % (non-condensing)				
Mains voltage	90 – 230 VAC, single phase, 47 – 63 Hz ¹⁴⁾				
Average power consumption	15 W	25 W	30	W	

- Due to continuous improvements all specifications are subject to change. Unless stated otherwise all specifications are measured at fundamental wavelength and maximum pulse repetition rate. The parameters marked typical are not specifications. They are indications of typical performance and will vary with each unit we manufacture.
- 2) Factory-set pulse repetition rate is fixed at max repetition rate shown in the table. Smoothly variable pulse repetition rate is possible when laser is externally triggered. In internal triggering mode repetition rate can be divided by integer number down to f/2, f/3, f/4,... 1 Hz.
- 3) At FWHM level at fundamental wavelength, measured with 350 ps rise time photodiode.
- When pul se energy is presented in xx/yy format, first number is for 1064 nm wavelength version, second – for 1053 nm version. Pulse energies for harmonic wavelengths are maximum possible values available from H1 series harmonic generator module. See datasheets of H1 for details.
- Measured during 30 seconds operation after warm-up.
- 6) Over 8 hour period after 20 minutes of warm-up when ambient temperature variation is less than ±2 °C.
- ⁷⁾ Full angle measured at the 4σ level.
- Peak to peak value, measured for 30 seconds of operation after laser warm-up.
- 9) Beam diameter is measured 20 cm from laser output at the 4σ level.
- ¹⁰⁾ In respect to falling edge of pump diode triggering pulse.
- 11) Motorized attenuator intended to be attached to the laser housing. Transmission can be changed remotely trough laser web-server control interface.
- Power adapter dimensions might differ from indicated here, depending on model.
- 13) Q1 laser is compatible all models of stand-alone H1 harmonics generator. Please refer to harmonic generator datasheets for detailed specifications.
- 14) Laser can be powered from appropriate 12 VDC power source. Please inquire for details.








SPECIFICATIONS AT 20 Hz PULSE REPETITION RATE 1)

MODEL.	Q1				
MODEL	-B20	-C20	-D20		
Wavelength	1064 or 1053 nm				
Pulse repetition rate 2)	single-shot to 20 Hz				
Pulse duration 3)	< 8 ns				
Pulse energy 4)					
1064 / 1053 nm	10 mJ	15 mJ	25 mJ		
532 / 526.5 nm	4 mJ	7.5 mJ	12.5 mJ		
355 / 351 nm	2.5 mJ	5 mJ	7.5 mJ		
266 / 263 nm	1 mJ	1.8 mJ	3 mJ		
213 / 211 nm	0.4 mJ	0.7 mJ	1.2 mJ		
Pulse to pulse energy stability 5)					
1064 / 1053 nm	< 0.5 % RMS				
532 / 526.5 nm	< 2.5 % RMS				
355 / 351 nm	< 3.5 % RMS				
266 / 263 nm	< 4 % RMS				
213 / 211 nm	< 5 % RMS				
Power drift 6)	± 3.0 %				
Beam profile	bell-shaped, >75 % fit to Gaussian				
Beam divergence 7)	< 1.5 mrad < 1 mrad				
Beam pointing stability 8)	< 20 μrad				
Polarization		linear, horizontal			
Typical beam diameter 9)	1.7 mm	2.0 mm	2.5 mm		
Jitter 10)	< 0.5 ns RMS				
OPTIONAL ATTENUATOR 11)					
Transmission range	1 – 95 %				
Transmission range		1 - 73 /0			
DIMENSIONS					
Laser head (W×L×H)	$110 \times 231 \times 112 \text{ mm}^3$				
Controller unit (W×L×H)	$108 \times 191 \times 59 \text{ mm}^3$				
Power adapter (W×L×H) 12)	50 × 125 × 31 mm ³				
Harmonic generator (W×L×H) 13)	$110 \times 242 \times 112 \text{ mm}^3$				
OPERATING REQUIREMENTS					
Cooling requirements	air cooled				
Ambient temperature	15 – 30 °C				
Relative humidity	10 – 80 % (non-condensing)				
Mains voltage	90 – 230 VAC, single phase, 47 – 63 Hz ¹⁴⁾				
Average power consumption	30 W				

- Due to continuous improvements all specifications are subject to change. Unless stated otherwise all specifications are measured at fundamental wavelength and maximum pulse repetition rate. The parameters marked typical are not specifications. They are indications of typical performance and will vary with each unit we manufacture.
- ²⁾ Factory-set pulse repetition rate is fixed at max repetition rate shown in the table. Smoothly variable pulse repetition rate is possible when laser is externally triggered. In internal triggering mode repetition rate can be divided by integer number down to f/2, f/3, f/4,... 1 Hz.
- 3) At FWHM level at fundamental wavelength, measured with 350 ps rise time photodiode.
- When pul se energy is presented in xx/yy format, first number is for 1064 nm wavelength version, second – for 1053 nm version. Pulse energies for harmonic wavelengths are maximum possible values available from H1 series harmonic generator module. See datasheets of H1 for details.
- Measured during 30 seconds operation after warm-up.
- 6) Over 8 hour period after 20 minutes of warm-up when ambient temperature variation is less than ±2 °C.
- ⁷⁾ Full angle measured at the 4σ level.
- Peak to peak value, measured for 30 seconds of operation after laser warm-up.
- 9) Beam diameter is measured 20 cm from laser output at the 4σ level.
- 10) In respect to falling edge of pump diode triggering pulse.
- 11) Motorized attenuator intended to be attached to the laser housing. Transmission can be changed remotely trough laser web-server control interface.
- 12) Power adapter dimensions might differ from indicated here, depending on model.
- 13) Q1 laser is compatible all models of stand-alone H1 harmonics generator. Please refer to harmonic generator datasheets for detailed specifications.
- 14) Laser can be powered from appropriate 12 VDC power source. Please inquire for details.











Laser head Q1 with controller unit



SPECIFICATIONS AT 50 Hz PULSE REPETITION RATE 1)

MODEL	Q1				
MODEL	-50	-A50	-B50		
Wavelength	1064 or 1053 nm				
Pulse repetition rate ²⁾	single-shot to 50 Hz				
Pulse duration 3)	< 10 ns	< 9 ns	< 8 ns		
Pulse energy 4)					
1064 / 1053 nm	2.5 mJ	5 mJ	10 mJ		
532 / 526.5 nm	1.2 mJ	2.5 mJ	5 mJ		
355 / 351 nm	0.6 mJ	1.5 mJ	3 mJ		
266 / 263 nm	0.3 mJ	0.7 mJ	1.5 mJ		
213 / 211 nm	0.1 mJ	0.2 mJ	0.5 mJ		
Pulse to pulse energy stability 5)					
1064 / 1053 nm	< 0.5 % RMS				
532 / 526.5 nm	< 2.5 % RMS				
355 / 351 nm	< 3.5 % RMS				
266 / 263 nm	< 4 % RMS				
213 / 211 nm	< 5 % RMS				
Power drift ⁶⁾	± 3.0 %				
Beam profile	nearly TEM ₀₀ , >80 % fit to Gaussian				
Beam divergence 7)	< 1.5 mrad				
Beam pointing stability 8)	< 20 μrad				
Polarization	linear, horizontal				
Typical beam diameter 9)	1.5 mm				
Jitter 10)	< 0.5 ns RMS				
OPTIONAL ATTENUATOR 11)					
Transmission range	1 – 95 %				
DIMENSIONS					
Laser head (W×L×H)	110 × 231× 112 mm ³				
Controller unit (W×L×H)	108 × 191 × 59 mm ³				
Power adapter (W×L×H) 12)	$50 \times 125 \times 31 \text{ mm}^3$				
Harmonic generator (W×L×H) 13)	110 × 242 × 112 mm³				
OPERATING REQUIREMENTS					
Cooling requirements	air cooled				
Ambient temperature	15 – 30 °C				
Relative humidity	10 – 80 % (non-condensing)				
Mains voltage	90 – 230 VAC, single phase, 47 – 63 Hz ¹⁴⁾				
Average power consumption	<40 W				

- Due to continuous improvements all specifications are subject to change. Unless stated otherwise all specifications are measured at fundamental wavelength and maximum pulse repetition rate. The parameters marked typical are not specifications. They are indications of typical performance and will vary with each unit we manufacture.
- ²⁾ Factory-set pulse repetition rate is fixed at max repetition rate shown in the table. Smoothly variable pulse repetition rate is possible when laser is externally triggered. In internal triggering mode repetition rate can be divided by integer number down to f/2, f/3, f/4,.. 1 Hz.
- ³⁾ At FWHM level at fundamental wavelength, measured with 350 ps rise time photodiode.
- When pul se energy is presented in xx/yy format, first number is for 1064 nm wavelength version, second – for 1053 nm version. Pulse energies for harmonic wavelengths are maximum possible values available from H1 series harmonic generator module. See datasheets of H1 for details.
- Measured during 30 seconds operation after warm-up.
- 6) Over 8 hour period after 20 minutes of warm-up when ambient temperature variation is less than ±2 °C.
- $^{7)}$ Full angle measured at the 4σ level.
- Peak to peak value, measured for 30 seconds of operation after laser warm-up.
- 9) Beam diameter is measured 20 cm from laser output at the 4σ level.
- ¹⁰⁾ In respect to falling edge of pump diode triggering pulse.
- 11) Motorized attenuator intended to be attached to the laser housing. Transmission can be changed remotely trough laser web-server control interface.
- Power adapter dimensions might differ from indicated here, depending on model.
- 13) Q1 laser is compatible all models of stand-alone H1 harmonics generator. Please refer to harmonic generator datasheets for detailed specifications.
- 14) Laser can be powered from appropriate 12 VDC power source. Please inquire for details.

AUXILIARY EQUIPMENT

Stand-alone H1 series harmonic generators for down to 211 nm output

Attachable motorized attenuator for fundamental wavelength beam

Attachable pulse energy monitor with analog and/or digital output


Baseplate for mounting Q1 and H1 or Q1 and fiber coupler together


Stand-alone air-purging unit for long lifetime of UV optics











Laser head with attachable motorized attenuator dimensions (in mm)



Quantum Light Instruments Ltd. Mokslininku 6A LT-08412, Vilnius, Lithuania Phone: +370 5 250 3717 Fax: +370 5 250 3716

Email: sales@qlinstruments.com