



High Etendue Spectrometer

Technical Specifications

Contact Details

Correspondence
Address:

IS-Instruments Ltd
Pipers Business Centre
220 Vale Road
Tonbridge
Kent
TN9 1SP
UK

Email:

info@is-instruments.com

Phone:

+44 (0) 1732 373 020

Fax:

+44 (0) 1732 373 001



Technical Specifications

This document shows the technical specifications of the HES range. Further information can be found on the IS-Instruments website, where you can download gold standard scientific articles of the instrument performance when used in the field. The HES range is predominately used for Raman spectroscopy, and the laser and probe are also provided for a complete Raman unit. Bespoke systems are available on request.

Model	HES1000	HES2000	HES2003	HES2000IR
Configuration	SHS	SHS	SHS	SHS (Infra-Red)
Wavelength Range	350 – 1000 nm units available Typical setup for Raman (@532)	350 – 1000 nm units available Typical setup for Raman (@ 785 nm)	350 – 1000 nm units available Typical setup for Raman (@785 nm)	850-2000 nm unit available Typical setup for Raman (@1064 nm)
Range can be adjusted as required from 50 – 4000 cm ⁻¹	200-2500cm ⁻¹	200-2500cm ⁻¹	200-2500cm ⁻¹	200-2500cm ⁻¹
Resolution (per Fourier bin)	< 3cm ⁻¹			< 6cm ⁻¹
< 1 cm ⁻¹ (also available)				
Slit	No Slit			
Fibre Input	SMA FC/PC		Custom	SMA FC/PC
Fibre diameter	1 mm		3 mm	1 mm (3 mm also available)
Fibre NA	0.22			
Linearity	> 99 %			
Detector Type	Machine vision camera (max integration time 10 seconds)	Cooled CCD	Cooled CCD	Cooled InGaAs CCD
Supply Voltage				
Dimensions				
Weight				
Software				

Example spectra

The most common configuration for the HES instrument is as the main workhorse within a Raman Spectrometer. The Spatial Heterodyne Spectrometer (SHS) configuration provides the system with greater than x100 throughput, allowing the system to be used for standoff Raman measurements, as well as for Raman observations of diffuse targets. Below are examples of Raman spectra observed with a HES2000 instrument.




Figure 1: Raman spectra: (a) Aspirin and (b) Paracetamol

The HES range of instruments can be constructed with any detector, however ISI's detector of choice is a Cooled CCD from Andor. The Quantum efficiency of these detectors as a function of wavelength is given below.








Figure 2: Quantum efficiency of the cooled CCD

Mechanical Drawings



Disclaimer

THE INFORMATION CONTAINED HEREIN IS PROVIDED "AS IS" WITHOUT WARRANTY, CONDITION OR REPRESENTATION OF ANY KIND, EITHER EXPRESS, IMPLIED, STATUTORY OR OTHERWISE, INCLUDING BUT NOT LIMITED TO, ANY WARRANTY OF MERCHANTABILITY, NON-INFRINGEMENT OR FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT SHALL IS-INSTRUMENTS BE LIABLE FOR ANY LOSS OR DAMAGE, WHETHER DIRECT, INDIRECT, SPECIAL, INCIDENTAL, CONSEQUENTIAL OR OTHERWISE HOWSOEVER CAUSED WHETHER ARISING IN CONTRACT TORT OR OTHERWISE, ARISING OUT OF OR IN CONNECTION WITH THE USE OF THE INFORMATION PROVIDED HEREIN.

COPYRIGHT AND PROTECTIVE NOTICES:

The copyright in this document and the associated drawings are the property of IS-Instruments Ltd. and all rights are reserved. This document and the associated drawings are issued on condition that they are not copied, reprinted, or reproduced, nor their contents disclosed. The publication of information in this documentation does not imply freedom from any patent or proprietary right of IS-Instruments or any third party.

TRADEMARKS & PATENT INFORMATION:

IS-Instruments and the ISI logo are trademarks of IS-instruments Ltd. All other marks are property of their owners