


FTSS 355-50

Diode pumped passively Q-switched solid state laser

- 355 nm
- single pulse
- < 1.3 ns
- 1 – 100 Hz
- > 70 µJ




biology · biomedicine · chemistry · analytics

Optical Data	Wavelength	355 nm
	Spatial Mode	TEM ₀₀
	M ²	< 1.5
	Beam Divergence (full angle)	< 3.5 mrad
	Beam Ellipticity	< 1.5 : 1
	Waist Diameter	280 ± 80 µm (located at about 80mm inside the laser head)
	Beam Diameter	300 ± 100 µm (at laser exit)
	Peak Power	> 60 kW @ 1 - 100 Hz
	Pulse Energy	> 70 µJ @ 1 - 100 Hz
	Pulse Repetition Rate (with external trigger)	1 - 100 Hz
	Pulse Width (FWHM)	< 1.3 ns
	Polarization Ratio	> 100:1, vertical
	Long term pulse energy stability (6 hours) ¹⁾	< ± 3 %
	Pulse-To-Pulse Stability ²⁾	< 1 % rms
	Laser Classification	3B / IIIb
	Residual Emission	< 0.5 mW @ 1064 nm (Class 3R) < 0.1 mW @ 532 nm (Class 2)
Optical Output	Free Beam	
Electrical Data	Electrical Power Consumption	< 70 W
	Line Voltage	100 - 240 V AC (50-60 Hz) or 24 V DC
Interface	RS 232, USB	
Miscellaneous	Warm-up Time	< 15 min
	Operating Temperature	18 - 38 °C
	Laser Head Size	217 x 65 x 45 mm (core dimensions)
Options	SMA-connector for fiber with core diameter ≥ 100 µm	
	Synchronization signal output (rise time < 2 ns)	
	Manual shutter or electrical beam blocker	
	Manual or electrical driven wavelength switch 355 nm / 532 nm	
	Higher repetition rates on request	
	External telescope (e.g. M=5)	
	Manual or electrical attenuator	
Stand Alone System (incl. key-switch, heat sink, manual shutter; CDRH compliant)		

¹⁾ Drift over 6 hours, energy averaged over 10 sec after 5 min of continuous operation, temperature variation < 3 °C/hour.
²⁾ RMS over 1000 pulses after 5 min of continuous operation.

Laser Head and Controller



Laser Safety Labels

The FTSS355-50 lasers are class 3B / IIIb according to IEC 60825-1:2014



© Oct. 2021, Rev 2.5, CryLaS GmbH. All rights reserved.
Product specifications are subject to change without notice.

