

Precise Focusing Stage



Accessory can be combined with all models of miniRaman spectrometer

- Provides protection from ambient light
- Provides a safe operation since user will not see the laser irradiation
- Provides precise focusing of the laser beam on a sample

Precise Focusing Stage XYZ



- lateral travel range: 25 mm along the X and Y axes
- axial travel range: 25 mm
- precise focusing of the laser beam on a sample in axial and lateral dimensions with accuracy 10 μm
- magnetic fixation of the stage cover in up and down positions
- weight: 3,2 kg
- dimensions: 215 mm x 175 mm x 320,5 mm

Precise Focusing Stage Z



- axial travel range: 25 mm
- precise focusing of the laser beam on a sample in axial dimension with accuracy 10 μm
- magnetic fixation of the stage cover in up and down positions
- weight: 2,5 kg
- dimensions: 215 mm x 175 mm x 320,5 mm



Lightnovol ApS
Blokken 15, 1. tv.
3460 Birkerød
Denmark (DK)
CVR: 40979603

+45 71 37 04 10
info@lightnovol.com
<https://lightnovol.com>

